

## Danaidal

## Chemical Properties


CAS No. :

Formula: C<sub>8</sub>H<sub>9</sub>NO

Molecular Weight: 135.166

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Danaidal is a useful organic compound for research related to life sciences and the catalog number is T125138.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.3981 mL	36.9905 mL	73.9809 mL
5 mM	1.4796 mL	7.3981 mL	14.7962 mL
10 mM	0.7398 mL	3.699 mL	7.3981 mL
50 mM	0.148 mL	0.7398 mL	1.4796 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481