


2',3',5'-Tri-O-benzoyl-5-iodouridine

Chemical Properties

CAS No. : 2880-91-3
 Formula: C₃₀H₂₃IN₂O₉
 Molecular Weight: 682.42
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nucleoside Derivatives - Halo-Nucleosides, 5-Modified pyrimidine nucleosides
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4654 mL	7.3269 mL	14.6537 mL
5 mM	0.2931 mL	1.4654 mL	2.9307 mL
10 mM	0.1465 mL	0.7327 mL	1.4654 mL
50 mM	0.0293 mL	0.1465 mL	0.2931 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481