


1H-pyrazolo[3,4-b]pyridine,3-methyl-

Chemical Properties

CAS No. :	116834-96-9
Formula:	C7H7N3
Molecular Weight:	133.154
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1H-pyrazolo[3,4-b]pyridine,3-methyl-, with CAS No. 116834-96-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1H-pyrazolo[3,4-b]pyridine,3-methyl- provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.5103 mL	37.5516 mL	75.1033 mL
5 mM	1.5021 mL	7.5103 mL	15.0207 mL
10 mM	0.751 mL	3.7552 mL	7.5103 mL
50 mM	0.1502 mL	0.751 mL	1.5021 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481