


Madrasin

Chemical Properties

CAS No. :	374913-63-0
Formula:	C ₁₆ H ₁₇ N ₅ O ₂
Molecular Weight:	311.34
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Madrasin (DDD00107587) is a potent and cell penetrant splicing inhibitor that interferes with the early stages of spliceosome assembly.
Targets(IC50)	Others,DNA/RNA Synthesis

Solubility Information

Solubility	DMSO: 1 mg/mL (3.21 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2119 mL	16.0596 mL	32.1192 mL
5 mM	0.6424 mL	3.2119 mL	6.4238 mL
10 mM	0.3212 mL	1.606 mL	3.2119 mL
50 mM	0.0642 mL	0.3212 mL	0.6424 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Andrea Pawellek et al. J Biol Chem, 2014 Dec 12, 289(50): 34683-34698.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481