

Compound T129970

Chemical Properties


CAS No. :

Formula: C₂₃H₃₇N₀O₈

Molecular Weight: 455.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-hexyl-3-methoxy-2-methyl-8-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-1,4,5,6,7,8-hexahydroquinolin-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1951 mL	10.9757 mL	21.9515 mL
5 mM	0.439 mL	2.1951 mL	4.3903 mL
10 mM	0.2195 mL	1.0976 mL	2.1951 mL
50 mM	0.0439 mL	0.2195 mL	0.439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481