

## Fmoc-Gly-Ser(psi(Me,Me)pro)-OH acetate

### Chemical Properties

CAS No. :


Formula: C<sub>25</sub>H<sub>28</sub>N<sub>2</sub>O<sub>8</sub>

Molecular Weight: 484.5

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	Fmoc-Gly-Ser(psi(Me,Me)pro)-OH acetate is a dipeptide.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.064 mL	10.3199 mL	20.6398 mL
5 mM	0.4128 mL	2.064 mL	4.128 mL
10 mM	0.2064 mL	1.032 mL	2.064 mL
50 mM	0.0413 mL	0.2064 mL	0.4128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481