

Compound N004-0003

Chemical Properties


CAS No. :

Formula: C₃₂H₄₅N₂O₆

Molecular Weight: 539.713

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N004-0003 is a useful organic compound for research related to life sciences and the catalog number is T131381.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8528 mL	9.2642 mL	18.5285 mL
5 mM	0.3706 mL	1.8528 mL	3.7057 mL
10 mM	0.1853 mL	0.9264 mL	1.8528 mL
50 mM	0.0371 mL	0.1853 mL	0.3706 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481