

Compound N039-0038

## Chemical Properties


CAS No. :

Formula: C45H74O17

Molecular Weight: 887.07

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N039-0038 is a useful organic compound for research related to life sciences and the catalog number is T131267.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1273 mL	5.6365 mL	11.2731 mL
5 mM	0.2255 mL	1.1273 mL	2.2546 mL
10 mM	0.1127 mL	0.5637 mL	1.1273 mL
50 mM	0.0225 mL	0.1127 mL	0.2255 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481