

Diethyl malonate

Chemical Properties


CAS No. : 105-53-3

Formula: C7H12O4

Molecular Weight: 160.17

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Diethyl malonate is naturally occurring in grapes and strawberries.
Targets(IC50)	Others

Solubility Information

Solubility	H2O: 10 mg/mL (62.43 mM), Sonication is recommended. DMSO: 60 mg/mL (374.6 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2434 mL	31.2168 mL	62.4337 mL
5 mM	1.2487 mL	6.2434 mL	12.4867 mL
10 mM	0.6243 mL	3.1217 mL	6.2434 mL
50 mM	0.1249 mL	0.6243 mL	1.2487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481