

3-Methoxyshancigusin I

Chemical Properties

CAS No. :

Formula: C₂₉H₃₆O₁₅

Molecular Weight: 624.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-Methoxyshancigusin I can be used in related research in the field of life sciences, and its product number is TN6576.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.601 mL	8.0051 mL	16.0102 mL
5 mM	0.3202 mL	1.601 mL	3.202 mL
10 mM	0.1601 mL	0.8005 mL	1.601 mL
50 mM	0.032 mL	0.1601 mL	0.3202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481