


Ferulamide

Chemical Properties

CAS No. :	61012-31-5
Formula:	C ₁₀ H ₁₁ NO ₃
Molecular Weight:	193.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ferulamide (4-Hydroxy-3-methoxycinnamide), a Ferulic acid derivative, shows potent inhibitory activity against arachidonic acid-induced platelet aggregation.
Targets(IC50)	AChR,PAFR
In vitro	A series of Ferulamide derivatives were prepared and evaluated for their anti-platelet activities. Some of these compounds showed potent inhibitory activity against arachidonic acid-induced platelet aggregation. Their structure-activity relationships are also discussed[1].

Solubility Information

Solubility	DMSO: 50 mg/mL (258.8 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.176 mL	25.8799 mL	51.7598 mL
5 mM	1.0352 mL	5.176 mL	10.352 mL
10 mM	0.5176 mL	2.588 mL	5.176 mL
50 mM	0.1035 mL	0.5176 mL	1.0352 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Synthesis and anti-platelet activity of ferulamide derivatives *Journal of Chinese Pharmaceutical Sciences*, 2002, 54 (1):41-52.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481