


## Ylangenol

## Chemical Properties

CAS No. :	41610-69-9
Formula:	C <sub>15</sub> H <sub>24</sub> O
Molecular Weight:	220.35
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ylangenol is a natural product of <i>Laggera</i> , Asteraceae. The catalog number is TN5361 and the CAS number is 41610-69-9. Ylangenol can be used as a reference standard.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5382 mL	22.6912 mL	45.3823 mL
5 mM	0.9076 mL	4.5382 mL	9.0765 mL
10 mM	0.4538 mL	2.2691 mL	4.5382 mL
50 mM	0.0908 mL	0.4538 mL	0.9076 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481