


2-Phenylbutanoic acid

Chemical Properties

CAS No. :	90-27-7
Formula:	C ₁₀ H ₁₂ O ₂
Molecular Weight:	164.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Phenylbutanoic acid (2-Phenylbutyric acid) is used as an industrial intermediate.
Targets(IC50)	HIV Protease,Endogenous Metabolite

Solubility Information

Solubility	DMSO: 10 mg/mL (60.9 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (6.09 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0901 mL	30.4507 mL	60.9013 mL
5 mM	1.218 mL	6.0901 mL	12.1803 mL
10 mM	0.609 mL	3.0451 mL	6.0901 mL
50 mM	0.1218 mL	0.609 mL	1.218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

MATHIVAT A, COTTET J. Bull Mem Soc Med Hop Paris. 1953 Nov 27-Dec 4;69(30-31):1030-48.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481