


Cryptomerin B

Chemical Properties

CAS No. : 22012-98-2
 Formula: C₃₂H₂₂O₁₀
 Molecular Weight: 566.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cryptomerin B is a natural product for research related to life sciences. The catalog number is TN5510 and the CAS number is 22012-98-2.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7652 mL	8.826 mL	17.6519 mL
5 mM	0.353 mL	1.7652 mL	3.5304 mL
10 mM	0.1765 mL	0.8826 mL	1.7652 mL
50 mM	0.0353 mL	0.1765 mL	0.353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481