


## 5,5'-Dimethoxylariciresinol

## Chemical Properties

CAS No. :	116498-58-9
Formula:	C <sub>22</sub> H <sub>28</sub> O <sub>8</sub>
Molecular Weight:	420.458
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5,5'-Dimethoxylariciresinol is a natural product of Phellodendron, Rutaceae. The catalog number is TN3078 and the CAS number is 116498-58-9. 5,5'-Dimethoxylariciresinol can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3783 mL	11.8917 mL	23.7835 mL
5 mM	0.4757 mL	2.3783 mL	4.7567 mL
10 mM	0.2378 mL	1.1892 mL	2.3783 mL
50 mM	0.0476 mL	0.2378 mL	0.4757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481