

Acid-PEG2-C2-Boc

Chemical Properties

CAS No. : 2086688-99-3


Formula: C₁₂H₂₂O₆

Molecular Weight: 262.3

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acid-PEG2-C2-Boc is a PEG-based linker for PROTACs, joining two essential ligands crucial for forming PROTAC molecules, enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two different ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, and the other targets the protein of interest. These compounds utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8124 mL	19.0621 mL	38.1243 mL
5 mM	0.7625 mL	3.8124 mL	7.6249 mL
10 mM	0.3812 mL	1.9062 mL	3.8124 mL
50 mM	0.0762 mL	0.3812 mL	0.7625 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481