


## cis-Martynoside

## Chemical Properties

CAS No. :	155899-92-6
Formula:	C <sub>31</sub> H <sub>40</sub> O <sub>15</sub>
Molecular Weight:	652.64
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	cis-Martynoside is a natural product for research related to life sciences. The catalog number is TN5745 and the CAS number is 155899-92-6.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5322 mL	7.6612 mL	15.3224 mL
5 mM	0.3064 mL	1.5322 mL	3.0645 mL
10 mM	0.1532 mL	0.7661 mL	1.5322 mL
50 mM	0.0306 mL	0.1532 mL	0.3064 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481