

2-(2-acetyl-3,5-dihydroxy-4,6-dimethylphenyl)acetic acid

Chemical Properties


CAS No. :

Formula: C₁₂H₁₄O₅

Molecular Weight: 238.24

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

2-(2-acetyl-3,5-dihydroxy-4,6-dimethylphenyl)acetic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1974 mL	20.9872 mL	41.9745 mL
5 mM	0.8395 mL	4.1974 mL	8.3949 mL
10 mM	0.4197 mL	2.0987 mL	4.1974 mL
50 mM	0.0839 mL	0.4197 mL	0.8395 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481