

PIK-124

## Chemical Properties


CAS No. : 711025-68-2

Formula: C<sub>19</sub>H<sub>12</sub>ClFN<sub>2</sub>O<sub>3</sub>S<sub>2</sub>

Molecular Weight: 434.89

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PIK-124 is a potent PI3K inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2994 mL	11.4972 mL	22.9943 mL
5 mM	0.4599 mL	2.2994 mL	4.5989 mL
10 mM	0.2299 mL	1.1497 mL	2.2994 mL
50 mM	0.046 mL	0.2299 mL	0.4599 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481