

NTU281

## Chemical Properties


CAS No. : 815619-12-6

Formula: C<sub>25</sub>H<sub>31</sub>N<sub>2</sub>O<sub>6</sub>S

Molecular Weight: 487.59

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	NTU281 is a potent inhibitor of transglutaminase-2 that reduces the elevation of serum creatinine and proteinuria in diabetic rats and decreases glomerular type I collagen accumulation, Hic-5 and $\alpha$ -SMA expression, as well as apoptosis, and can be used in studies of diabetes-induced glomerulosclerosis.
Targets(IC50)	Apoptosis,Others,Glutaminase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0509 mL	10.2545 mL	20.509 mL
5 mM	0.4102 mL	2.0509 mL	4.1018 mL
10 mM	0.2051 mL	1.0255 mL	2.0509 mL
50 mM	0.041 mL	0.2051 mL	0.4102 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481