


## 5'-O-Triphenylmethyl-2'-deoxyuridine

### Chemical Properties

CAS No. : 14270-73-6  
 Formula: C<sub>28</sub>H<sub>26</sub>N<sub>2</sub>O<sub>5</sub>  
 Molecular Weight: 470.52  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	5'-O-Triphenylmethyl-2'-deoxyuridine is a Nucleoside Derivative - Protected nucleoside with NH <sub>2</sub> /OH open.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1253 mL	10.6265 mL	21.2531 mL
5 mM	0.4251 mL	2.1253 mL	4.2506 mL
10 mM	0.2125 mL	1.0627 mL	2.1253 mL
50 mM	0.0425 mL	0.2125 mL	0.4251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481