

## AG-3 Protein, Human, Recombinant (His)

### General Information

Synonyms:	BCMP11;hAG-3;AG-3;AG3;Anterior Gradient Protein 3 Homolog;Breast Cancer Membrane Protein 11;AGR3
Protein Construction:	Ile22-Leu166
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q8TD06
Molecular Weight:	16 KDa (reducing condition)
AA Sequence:	Ile22-Leu166

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM Tris-HCl, 150 mM NaCl, 2 mM EDTA, pH 8.5.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Anterior Gradient Protein 2(AG-2) and Anterior Gradient Protein 3 (AG-3) are human homologues of genes involved in differentiation, are associated with oestrogen receptor-positive breast tumours and interact with metastasis gene C4.4a and dystroglycan (hAG-3 protein). AG-3 could serve as a prognostic marker for survival in patients with low grade and high grade serous ovarian carcinomas.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481