

Galectin-9 Protein, Human, Recombinant (GST)

General Information

Synonyms:	Gal-9delta5;Ecalectin;LGALS9;Galectin-9;Gal-9δ5;Gal-9;Tumor antigen HOM-HD-21
Protein Construction:	Met1-Thr323
Species:	Human
Expression Host:	E. coli
Accession:	O00182-2
Molecular Weight:	60 KDa (reducing condition)
AA Sequence:	Met1-Thr323

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 1M PBS, 100 mM GSH, 15% Glycerol, pH 8.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Galectin-9 is a cytoplasmic protein that contains two galectin domains. Galectin-9 is an S-type lectin that is over-expressed in Hodgkin's disease tissue. Galectin-9 binds galactosides and has high affinity for the Forssman pentasaccharide. Galectin-9 plays a role in thymocyte-epithelial interactions relevant to the biology of the thymus and inhibits cell proliferation. Galectin-9 is a ligand for HAVCR2/TIM3 and induces T-helper type 1 lymphocyte (Th1) death. In addition, Galectin-9 suppresses tumor cell metastasis by interfering with the associations CD44, VCAM-1, Integrin α4β1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481