

Latent TGF-beta 2 Protein, Human, Recombinant (His & Avi), Biotinylated

General Information

Synonyms:	TGFB2;TGF- β -2;TGF-beta-2;Transforming growth factor beta-2 proprotein;Transforming growth factor β -2;Latent TGF- β 2;TGF β 2
Protein Construction:	Leu30-Ser390 (Cys33Ser)
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P61812
Molecular Weight:	13&39-50 KDa (reducing condition)
AA Sequence:	Leu30-Ser390

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Transforming growth factor beta (TGF β) is a multifunctional cytokine that regulates cell growth, differentiation, adhesion, migration and death dependent on cell type, developmental stage, or tissue conditions. There are three isoforms of TGF β (TGF β -1, -2 and -3). TGFB2 could be a potential biomarker for screening, surveillance and prognosis in GC, and that TGFB2 might drive EMT in GC through the NGF/TrKs pathway. Overexpression of TGFB2 may drive EMT through the neuriterelated signaling and affect TMB levels by regulating the DNA damage repair

A DRUG SCREENING EXPERT

pathways and immune infiltrates. TGF- β 2 is an immune suppressor involved in the development of immune tolerance, and recombinant TGF- β 2 incubation is more potent than TGF- β 1 or TGF- β 3 incubation in suppressing macrophage inflammatory responses.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481