

Fumarase Protein, Human, Recombinant

General Information

| | |
|-----------------------|--|
| Synonyms: | FH;Fumarase;Fumarate Hydratase Mitochondrial |
| Protein Construction: | Ala45-Lys510 |
| Species: | Human |
| Expression Host: | E. coli |
| Accession: | P07954 |
| Molecular Weight: | 45 KDa (reducing condition) |
| AA Sequence: | Ala45-Lys510 |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first. |
| Purity: | Greater than 95% as determined by reducing SDS-PAGE. (QC verified) |
| Endotoxin: | < 0.1 ng/μg (1 EU/μg) as determined by LAL test. |
| Formulation: | Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, pH 8.0. |

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Fumarase is an enzyme that catalyze the reversible hydration/dehydration of fumarate to S-malate and is involved in the tricarboxylic acid or Krebs cycle. Fumarase exists in both form, cytosolic form and N-terminal extend mitochondrial form. The N-terminal extended form is targeted to the mitochondrion, where the removal of the extension is the same form as in the cytoplasm. Fumarase is thought to act as a tumor suppressor, which deficiency can lead to progressive encephalopathy, cerebral atrophy and development delay.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481