

Otolin-1 Protein, Human, Recombinant (His)

General Information

Synonyms:	Otolin-1;otolin 1;C1qTNF15;OTOL1
Protein Construction:	Lys24-Pro477
Species:	Human
Expression Host:	HEK293 Cells
Accession:	A6NHN0
Molecular Weight:	84-94 KDa (reducing condition)
AA Sequence:	Lys24-Pro477

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, 5% Trehalose, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Otolin (OTOL1), also known as C1qTNF15, is an approximately 65 kDa protein found in the otoconial membrane lining the cochlea and vestibular labyrinth of the inner ear. Collagen-like protein specifically expressed in the inner ear, which provides an organic scaffold for otoconia, a calcium carbonate structure in the saccule and utricle of the ear. It associates into multimers and disulfide-linked oligomers and also associates with other otoconial proteins including and Otoconin-90 (also known as PLA2G2A, PLA2L, and phospholipase A2 homolog) and

Cerebellin-1. It is extensively glycosylated and has multiple hydroxylated proline residues in the collagenous regions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481