

PPP1R1A Protein, Human, Recombinant (His)

General Information

Synonyms:	PPP1R1A;I-1;IPP-1;Protein Phosphatase Inhibitor 1;Protein Phosphatase 1 Regulatory Subunit 1A;IPP1
Protein Construction:	Met1-Val171
Species:	Human
Expression Host:	E. coli
Accession:	AAH22470.1
Molecular Weight:	18 KDa (reducing condition)
AA Sequence:	Met1-Val171

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 90% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl, 1 mM DTT, 50% Glycerol, pH 8.5.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Protein Phosphatase 1 Regulatory Subunit 1A (PPP1R1A) is an inhibitor of protein-phosphatase 1. PPP1R1A is a cellular regulator of eIF2 alpha phosphorylation. In hormonal control of glycogen metabolism, IPP-1 protein plays important function. Hormones can elevate intracellular cAMP level and elevate IPP-1 activity. PPP1R1A activation caused cAMP increase, cAMP control over proteins that are not directly phosphorylated by PKA following a rise in intracellular calcium. IPP-1 is inactivated by calcineurin (PP2B). Multiple domains in IPP-1 target cellular PP1 complexes.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481