

ZNHIT1 Protein, Human, Recombinant (His)

General Information

Synonyms:	ZNFN4A1;CGBP1;inc finger HIT domain-containing protein 1;cyclin-G1-binding protein 1;p18 hamlet;ZNHIT1;zinc finger protein subfamily 4A member 1
Protein Construction:	Met1-Val154
Species:	Human
Expression Host:	E. coli
Accession:	O43257
Molecular Weight:	21 KDa (reducing condition)
AA Sequence:	Met1-Val154

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of PBS, 50% Glycerol, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

ZNHIT1 belongs to the ZNHIT1 family and contains one HIT-type zinc finger. It can be phosphorylated on Thr by MAPK11 or MAPK14. ZNHIT1 is a component of the chromatin-remodeling SRCAP complex, which is composed of at least SRCAP, DMAP, RUVBL1, RUVBL2, ACTL6A, YEATS4, ACTR6 and ZNHIT1. ZNHIT1 may play a role in p53-mediated apoptosis induction.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481