

STX6 Protein, Human, Recombinant

General Information

Synonyms:	STX6;Syntaxin-6
Protein Construction:	Ser2-Gln234
Species:	Human
Expression Host:	E. coli
Accession:	O43752
Molecular Weight:	32 KDa (reducing condition)
AA Sequence:	Ser2-Gln234

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM PB, 150 mM NaCl, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Syntaxin-6 (STX6) is a single-pass type IV membrane protein, which belongs to the Syntaxin family. STX6 is mainly localized in the plasma membrane. STX6 contains one t-SNARE coiled-coil homology domain and involved in intracellular vesicle trafficking. When STX6 function is inhibited, internalization through caveolae is dramatically reduced, whereas other endocytic mechanisms are unaffected. It is reported that STX6 is necessary for proper expression of focal adhesion kinase and integrin alpha5 adhesion receptor.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481