

PDYN Protein, Human, Recombinant (His)

General Information

Synonyms:	PDYN;ADCA;Big Dyn;Dynorphin B;PENKB;SCA23
Protein Construction:	Asp21-Ala254
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P01213
Molecular Weight:	35 KDa (reducing condition)
AA Sequence:	Asp21-Ala254

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 4 mM HCL.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Proenkephalin-B(PDYN), belongs to the opioid neuropeptide precursor family. The N-terminal domain contains 6 conserved cysteines thought to be involved in disulfide bonding and/or processing. Leu-enkephalins, which is a type of Proenkephalin-B, compete with and mimic the effects of opiate drugs. They play a role in a number of physiologic functions, including pain perception and responses to stress. Dynorphin peptides differentially regulate the kappa opioid receptor. Dynorphin A has a typical opioid activity, it is 700 times more potent than Leu-

enkephalin.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481