

Granzyme B/GZMB Protein, Mouse, Recombinant (His)

General Information

Synonyms:	Gzmb;Cytotoxic cell protease 1;Fragmentin-2;CTLA-1;Granzyme B(G;H);Ctla1;CCP1
Protein Construction:	Gly19-Ser247
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	P04187
Molecular Weight:	35 KDa (reducing condition)
AA Sequence:	Gly19-Ser247

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of PBS, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Granzyme B(GZMB) contains 1 peptidase S1 domain and belongs to the peptidase S1 family. This enzyme is necessary for target cell lysis in cell-mediated immune responses. It cleaves after Asp and seems to be linked to an activation cascade of caspases (aspartate-specific cysteine proteases) responsible for apoptosis execution. The protein cleaves caspase-3, -7, -9 and 10 to give rise to active enzymes mediating apoptosis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481