

## Serpin E2 Protein, Human, Recombinant (His)

### General Information

Synonyms:	PN-1;PI-7;Protease nexin 1;PI7;GDN;Peptidase inhibitor 7;Serpin E2;SERPINE2;Protease nexin I;PN1;Glia-derived Nexin
Protein Construction:	Ser20-Pro397
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P07093-2
Molecular Weight:	50 KDa (reducing condition)
AA Sequence:	Ser20-Pro397

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl, pH 8.0.

### Preparation and Storage

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

#### Shipping:

Proteins are shipped with blue ice.

### Protein Background

Serpin E2 is a secreted protein that belongs to the serpin family. Serpin E2 is a serine protease inhibitor with activity toward thrombin, trypsin, and urokinase. Serpin E2 expression is weak or absent in all normal pancreas and chronic pancreatitis tissue. In contrast, it was strongly over-expressed in the majority of pancreatic carcinoma as well as gastric and colorectal cancer samples. Serpin E2 promotes neurite extension by inhibiting thrombin. It also can bind heparin. It has been shown that Serpin E2 is a novel target of ERK signaling involved in human colorectal tumorigenesis. It plays an important role in controlling male fertility because its knockout male mice show a marked impairment in fertility from the onset of sexual maturity and its abnormal expression is found in the semen of men with seminal dysfunction.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481