

## BPGM Protein, Human, Recombinant (His)

### General Information

Synonyms:	3-Diphosphoglycerate Mutase;2;DPGM;Bisphosphoglycerate Mutase;3-Bisphosphoglycerate Synthase;BPG-Dependent PGAM;3-Bisphosphoglycerate Mutase Erythrocyte;BPGM
Protein Construction:	Ser2-Lys259
Species:	Human
Expression Host:	E. coli
Accession:	P07738
Molecular Weight:	30 KDa (reducing condition)
AA Sequence:	Ser2-Lys259

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 1 mM DTT, pH 8.0.

### Preparation and Storage

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

Proteins are shipped with blue ice.

### Protein Background

Bisphosphoglycerate Mutase (BPGM) is a member of the Phosphoglycerate Mutase family and BPG-Dependent PGAM subfamily. BPGM is a multifunctional enzyme. BPGM catalyzes 2,3-DPG synthesis via its synthetase activity, and 2,3-DPG degradation via its phosphatase activity. It also has phosphoglycerate phosphomutase activity. BPGM plays a major role in regulating hemoglobin oxygen affinity by controlling the levels of 2,3-bisphosphoglycerate (2,3-BPG). Deficiency of BPGM increases the affinity of cells for oxygen and result in hemolytic anemia.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481