

NPPB Protein, Human, Recombinant (His)

General Information

Synonyms:	NPPB;γ-Brain Natriuretic Peptide;Gamma-Brain Natriuretic Peptide;Natriuretic Peptides B
Protein Construction:	His27-His134
Species:	Human
Expression Host:	E. coli
Accession:	P16860
Molecular Weight:	16 KDa (reducing condition)
AA Sequence:	His27-His134

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM PB, 150 mM NaCl, 1 mM EDTA, 1 mM DTT, 20% Glycerol, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Human Natriuretic peptides B acts as a cardiac hormone; it is associated with many biological actions, such as diuresis, natriuresis, vasorelaxation, which inhibits the secretion of rennin and aldosterone. It acts as a paracrine antifibrotic factor in the heart. Natriuretic peptides B can help restore the body balance of salt and water, improves the heart function. Natriuretic peptides B binds and stimulates the cGMP production of the NPR1 receptor and binds the clearance receptor NPR3.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481