

DTYMK Protein, Human, Recombinant (His)

General Information

Synonyms:	DTYMK;TMPK;dTMP kinase;Thymidylate kinase;CDC8;TYMK
Protein Construction:	Met1-Lys212
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P23919
Molecular Weight:	28 KDa (reducing condition)
AA Sequence:	Met1-Lys212

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM PB, 5% Sucrose, 0.05% Tween 80, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Thymidylate kinase (DTYMK) is a ubiquitous enzyme of about 25 kD which belongs to thymidylate kinase family. DTYMK is important in the dTTP synthesis pathway for DNA synthesis. It participated in the pyrimidine metabolism pathway and dTTP biosynthesis pathway. DTYMK catalyzes the conversion of dTMP to dTDP and catalyzes the phosphorylation of thymidine 5'-monophosphate (dTMP) to form thymidine 5'-diphosphate (dTDP) in the presence of ATP and magnesium. Structural and functional analyses suggest that the cDNA codes for authentic human dTMP kinase. The mRNA levels and enzyme activities corresponded to cell cycle progression and cell growth stages.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481