

## HSP40 Protein, Human, Recombinant (His)

### General Information

Synonyms:	HDJ1;HSPF1;DNAJ1;Heat Shock 40 kDa Protein 1;HSP40;hDj-1;DNAJB1;Heat Shock Protein 40; Dnaj Homolog Subfamily B Member 1;Dnaj Protein Homolog 1;Human Dnaj Protein 1
Protein Construction:	Gly2-Ile340
Species:	Human
Expression Host:	E. coli
Accession:	P25685
Molecular Weight:	38 KDa (reducing condition)
AA Sequence:	Gly2-Ile340

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, 1 mM EDTA, pH 7.4.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Dnaj Homolog Subfamily B Member 1 (DNAJB1) is a member of the heat shock protein family. Heat shock proteins (HSPs) are a highly conserved family of stress response proteins. HSPs function primarily as molecular chaperones, facilitating the folding of other cellular proteins, preventing protein aggregation, or targeting improperly folded proteins to specific degradative pathways. DNAJB1 regulates cellular processes by aiding in the

## A DRUG SCREENING EXPERT

---

folding, transport and assembly. DNAJB1 contains a J-domain which controls interaction with the ATPase domain of DnaK. DNAJB1 interacts with HSP70 and can stimulate its ATPase activity. In addition, DNAJB1 stimulates the association between HSC70 and HIP.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481