

## MTHFS Protein, Human, Recombinant (His)

## General Information

Synonyms:	5-formyltetrahydrofolate cyclo-ligase;10-methenyl-tetrahydrofolate synthetase;MTHFS;5;Methenyl-THF synthetase
Protein Construction:	Met1-Ala203
Species:	Human
Expression Host:	E. coli
Accession:	P49914
Molecular Weight:	28 KDa (reducing condition)
AA Sequence:	Met1-Ala203

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 200 mM NaCl, 1 mM DTT, 50% Glycerol, pH 8.0.

## Preparation and Storage

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

Proteins are shipped with blue ice.

## Protein Background

5-formyltetrahydrofolate cyclo-ligase (MTHFS) belongs to the 5-formyltetrahydrofolate cyclo-ligase family. It is an enzyme that catalyzes the conversion of 5-formyltetrahydrofolate to 5,10-methenyltetrahydrofolate, contributes to tetrahydrofolate metabolism. MTHFS helps regulate carbon flow through the folate-dependent one-carbon metabolic network that supplies carbon for the biosynthesis of purines, thymidine and amino acids. An increased activity of the encoded protein can result in an increased folate turnover rate and folate depletion.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481