

SGSH Protein, Human, Recombinant (His)

General Information

Synonyms:	SGSH;N-Sulphoglucosamine Sulphohydrolase;Sulfo-glucosamine Sulfamidase; Sulphamidase;HSS
Protein Construction:	Arg21-Leu502
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P51688
Molecular Weight:	63 KDa (reducing condition)
AA Sequence:	Arg21-Leu502

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl, 1 mM CaCl ₂ , 10% Glycerol, pH 7.5.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

N-Sulphoglucosamine Sulphohydrolase (SGSH) is an important member of the sulfatase family which is involved in the degradation of heparin sulfate. SGSH binds one calcium ion per subunit as a cofactor. SGSH catalyzes N-sulfo-D-glucosamine and H₂O to D-glucosamine and sulfate. SGSH deficiency is result in mucopolysaccharidosis type 3A (MPS3A), a recessive lysosomal storage disease characterized by neurological dysfunction but relatively mild somatic manifestations.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481