

RPS7 Protein, Human, Recombinant (His)

General Information

Synonyms:	RPS7;40S ribosomal protein S7
Protein Construction:	Met1-Leu194
Species:	Human
Expression Host:	E. coli
Accession:	P62081
Molecular Weight:	23 KDa (reducing condition)
AA Sequence:	Met1-Leu194

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Acetate, 8% Sucrose, 50 mM NaCl, 0.05% Tween 80, 0.05% Tween 20, pH 4.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

40S ribosomal protein S7(RPS7) belongs to the S7E family of ribosomal proteins. It is phosphorylated by NEK6 during post-translational modification. RPS7 is located in the cytoplasm, binds IPO9 with high affinity. It also can interact with NEK6. As is required for rRNA maturation and typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. The unnormal expression of RPS7 may cause Diamond-Blackfan anemia 8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481