

## ARG2 Protein, Human, Recombinant (His)

### General Information

Synonyms:	ARG2;Kidney-type arginase;Type II arginase;Arginase-2, mitochondrial;Non-hepatic arginase
Protein Construction:	His24-Gly330
Species:	Human
Expression Host:	E. coli
Accession:	P78540
Molecular Weight:	33 KDa (reducing condition)
AA Sequence:	His24-Gly330

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 50 mM HEPES, 150 mM NaCl, pH 7.5.

### Preparation and Storage

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

Proteins are shipped with blue ice.

### Protein Background

Arginase-2 (ARG2) is a member of the arginase family. Arginase is a manganese-containing enzyme which catalyzes the hydrolysis of arginine to ornithine and urea. ARG2 is highly expressed in kidney and prostate, not founded in the liver, heart and pancreas. ARG2 has been implicated in the regulation of the arginine/ornithine concentrations in the cell. ARG2 may take part in the regulation of extra-urea cycle arginine metabolism and in down-regulation of nitric oxide synthesis. The extrahepatic arginase functions to regulate L-arginine bioavailability to NO synthase.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481