

PAP Protein, Human, Recombinant (His)

General Information

Synonyms:	PAP1;PAP;PDGFA-Associated Protein 1;HASPP28;28 kDa Heat- and Acid-Stable Phosphoprotein;PDAP1;PDGF-Associated Protein
Protein Construction:	Met1-Lys181
Species:	Human
Expression Host:	E. coli
Accession:	Q13442
Molecular Weight:	30 KDa (reducing condition)
AA Sequence:	Met1-Lys181

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 100 mM NaCl, 0.1 mM PMSF, 2 mM DTT, pH 8.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Human PAP, also known as 28 kDa heat- and acid-stable phosphoprotein, PDGF-associated protein, PDGFA-associated protein 1, PDAP1, HASPP28, is a protein which belongs to the PDAP1 family. The encoded protein in rodents has been shown to bind PDGFA with a low affinity. PDGF-Associated Protein (PAP) is a phosphoprotein that may enhance PDGFA-stimulated cell growth in fibroblasts, but inhibits the mitogenic effect of PDGFB. PDAP1 expression is induced by TNF-alpha, and cells overexpressing PDAP1 show significantly less apoptosis on exposure to TNF-alpha.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481