

UCMA Protein, Human, Recombinant

General Information

Synonyms:	UCMA;Unique cartilage matrix-associated protein;Gla-rich protein;C10orf49
Protein Construction:	Ser65-Thr138
Species:	Human
Expression Host:	E. coli
Accession:	Q8WVF2
Molecular Weight:	14 KDa (reducing condition)
AA Sequence:	Ser65-Thr138

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of PBS, pH 7.2.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

C10orf49 is a secreted protein which encoded by the UCMA gene. It is a member of the UCMA family. C10orf49 is predominantly expressed in resting chondrocytes. It may be involved in the negative control of osteogenic differentiation of osteochondrogenic precursor cells in peripheral zones of fetal cartilage and at the cartilage-bone interface.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481