

SEMA3C Protein, Human, Recombinant (His)

General Information

Synonyms:	Semaphorin-E;SEMA3C;SEMAE;Semaphorin-3C;Sema E
Protein Construction:	Gly21-Ser738
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q99985
Molecular Weight:	75-120&27KD (reducing condition)
AA Sequence:	Gly21-Ser738

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM NaH ₂ PO ₄ , 150 mM NaCl, 0.1M Arginine, 0.1M Glu, 10% Glycerol, 0.01% Tween20, 5% Trehalose, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Semaphorin 3C is one of six Class 3 secreted semaphorins, which are potent chemorepellents, that play a role in axon and/or vascular guidance during development can be up-regulated in the process of tumor progression. Semaphorin 3C consists of a 20 amino acids signal sequence, a N-terminal Sema domain, a cysteine knot, a furine-type cleavage site, an Ig-like domain, and a C-terminal basic domain. Semaphorin 3C expressed highly in the heart, colon, skeletal muscle, small intestine, testis, ovary and prostate, expressed lowly in the other organs, including the brain.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481