

SULT4A1 Protein, Human, Recombinant

General Information

Synonyms:	hBR-STL;SULTX3;Nervous System Sulfotransferase;SULT4A1;ST4A1;Brain Sulfotransferase-Like Protein;hBR-STL-1;NST;Sulfotransferase 4A1
Protein Construction:	Met1-Leu284
Species:	Human
Expression Host:	E. coli
Accession:	Q9BR01
Molecular Weight:	34 KDa (reducing condition)
AA Sequence:	Met1-Leu284

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM PB, 150 mM NaCl, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Sulfotransferase 4A1 (ST4A1) is a member of the Sulfotransferase 1 family. ST4A1 is highly expressed in the cerebral cortex and frontal lobe, but no expression is detected in the pancreas. ST4A1 is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. ST4A1 acts on catecholamines and T4 in a manner that may not involve sulfonation. ST4A1 may have a role in the metabolism of drugs and neurotransmitters in the CNS. In addition, ST4A1 is related to schizophrenia.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481