

ACAT2 Protein, Human, Recombinant (His)

General Information

Synonyms:	ACTL;Acetyl-CoA transferase-like protein;Acetyl-CoA acetyltransferase cytosolic;ACAT2;Cytosolic acetoacetyl-CoA thiolase
Protein Construction:	Met1-Glu397
Species:	Human
Expression Host:	E. coli
Accession:	Q9BWD1
Molecular Weight:	60 KDa (reducing condition)
AA Sequence:	Met1-Glu397

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, 10% Glycerol, pH 8.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

ACAT2 is a cytoplasmic enzyme which belongs to the thiolase family. ACAT2 takes part in lipid metabolism, lipoprotein assembly, catalyzing cholesterol esterification in mammalian cells. It is responsible for the synthesis of cholesteryl esters which are part of lipoproteins containing apoB. ACAT2 deficiency contributes to severe mental retardation and hypotonus.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481