

PACSIN1 Protein, Human, Recombinant (His)

General Information

Synonyms:	Protein Kinase C and Casein Kinase Substrate in Neurons Protein 1;KIAA1379;PACSIN1
Protein Construction:	Met1-Ile444
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q9BY11
Molecular Weight:	50-67 KDa (reducing condition)
AA Sequence:	Met1-Ile444

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM PB, 150 mM NaCl, pH 7.4.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Protein Kinase C and Casein Kinase Substrate in Neurons Protein 1 (PACSIN1) belongs to the PACSIN family. PACSIN1 contains one FCH domain and one SH3 domain. PACSIN1 is highly expressed in the brain and at lower levels in the heart, pancreas, and liver. PACSIN1 may play a role in vesicle formation and transport. PACSIN1 has been shown to interact with DNM1, PACSIN3, Huntingtin, and PACSIN2. In addition, PACSIN1 is phosphorylated by casein kinase 2 (CK2) and protein kinase C (PKC).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481