

IL-21R Protein, Mouse, Recombinant (aa 20-236, His)

General Information

Synonyms:	Nilr;IL-21R;Interleukin-21 receptor;IL-21 receptor;IL21r
Protein Construction:	Cys20-Pro236
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	Q9JHX3
Molecular Weight:	34-38 KDa (reducing condition)
AA Sequence:	Cys20-Pro236

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Interleukin-21 receptor (IL-21R) is a type I transmembrane glycoprotein within the class I cytokine receptor family, type 4 subfamily. IL-21R is expressed mainly on B cells, NK cells, and activated T cells, but is also found on dendritic cells, alternatively activated macrophages, intestinal lamina propria fibroblasts and epithelial cells, and keratinocytes. Both IL-21 and IL-4 are necessary for efficient B cell IgG1 production and normal germinal center architecture. B cell IL-21 R engagement induces Blimp-1 (which mediates plasma cell differentiation), and is

important for memory responses. IL-21R engagement on mouse NK cells enhances their cytotoxic activity and IFN- γ production. IL-21R engagement on CD8+ T cells aids control of viral infection and tumor growth; IL-21R is also necessary for sufficient numbers of regulatory T cells to combat chronic inflammation. IL-21R expression is often up-regulated in allergic skin inflammation, systemic lupus erythematosus and diffuse large B cell lymphoma (DLBCL).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481