

XPNPEP3 Protein, Human, Recombinant (His)

General Information

Synonyms:	XPNPEP3;X-Pro Aminopeptidase 3;Aminopeptidase P3;APP3;Probable Xaa-Pro Aminopeptidase 3
Protein Construction:	Met1-Ser507
Species:	Human
Expression Host:	E. coli
Accession:	Q9NQH7
Molecular Weight:	65 KDa (reducing condition)
AA Sequence:	Met1-Ser507

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 25 mM Tris-HCL, 1 mM DTT, pH 7.3.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Probable Xaa-Pro Aminopeptidase 3 (XPNPEP3) is a member of the peptidase M24B family. XPNPEP3 has two isoforms and both are widely expressed. XPNPEP3 is localized in the Mitochondrion. XPNPEP3 catalyzes the release of any N-terminal amino acid, including proline, that is linked to proline, even from a dipeptide or tripeptide. Defects in XPNPEP3 are the cause of nephronophthisis-like nephropathy type 1 which is a disorder with features of nephronophthisis, a cystic kidney disease leading to end-stage renal failure.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481