

THSD1 Protein, Human, Recombinant (His)

General Information

Synonyms:	Transmembrane Molecule with Thrombospondin Module;TMTSP;Thrombospondin Type-1 Domain-Containing Protein 1;THSD1
Protein Construction:	Glu25-Ile361
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q9NS62-2
Molecular Weight:	84 KDa (reducing condition)
AA Sequence:	Glu25-Ile361

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Thrombospondin Type-1 Domain-Containing Protein 1 (THSD1) is a single-pass type I membrane protein. THSD1 contains a signal peptide and one TSP type-1 domain that is found in thrombospondin. THSD1 is a good novel candidate for TSG as it has been mapped to 13q14. Alternatively spliced transcript variants encoding distinct isoforms have been observed. THSD1 may be involved in the complement pathway.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481