

ATG3 Protein, Human, Recombinant

General Information

Synonyms:	ATG3;Ubiquitin-like-conjugating enzyme ATG3;APG3;Protein PC3-96;Autophagy-related protein 3;APG3-like;APG3L;hApg3
Protein Construction:	Met1-Met314
Species:	Human
Expression Host:	E. coli
Accession:	Q9NT62
Molecular Weight:	38 KDa (reducing condition)
AA Sequence:	Met1-Met314

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	Greater than 95% as determined by reducing SDS-PAGE. (QC verified)
Endotoxin:	< 0.1 ng/μg (1 EU/μg) as determined by LAL test.
Formulation:	Supplied as a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl, 10% Glycerol, pH 8.0.

Preparation and Storage

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Ubiquitin-like-conjugating enzyme ATG3 (ATG3), also known as Apg3L and Apg3p, functions as a regulatory component of autophagosome biogenesis necessary for autophagy. ATG3 exhibits 98% aa sequence identity with both its mouse and rat orthologs. It is widely expressed and has highly levels in heart, skeletal muscle, kidney, liver and placenta. As an E2-like enzyme, involves in autophagy and mitochondrial homeostasis. ATG3 catalyzes the conjugation of ATG8-like proteins to PE which is essential for autophagy. ATG3 also can catalyze the conjugation of ATG12 to itself which plays a role in mitochondrial homeostasis but not in autophagy.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481